




Hệ thống sắc ký lỏng khối phổ Model: LC/MS 2020 (tiếp theo)

Tính kết quả



Chọn thư mục chứa data file





Bấm chuột phải vào mạng hình

Data View Parameters

Fragment Table MIC Settings Spectrum Format Delay Time

Segment# 1 Event# 1

☐ IIC
☐ BPC
☐ MIC 1
☒ None

	Disp.	Event	Factor
1	☒	1	1.00

	Disp.	Event	m/z	Factor
1	☒	1	254.00	1.00
2	☐			0.00
3	☐			0.00
4	☐			0.00
5	☐			0.00
6	☐			0.00
7	☐			0.00
8	☐			0.00
9	☐			0.00
10	☐			0.00

☐ Set factor on registering MC

☒ Base Shift

OK Cancel Apply Help

Chọn số khối cần hiển thị

Data View Parameters

Fragment Table MIC Settings Spectrum Format Delay Time

Segment# 1 Event# 1

☐ IIC
☐ BPC
☐ MIC 1
☒ None

	Disp.	Event	Factor
1	☒	1	1.00

	Disp.	Event	m/z	Factor
1	☒	1	254.00	1.00
2	☐			0.00
3	☐			0.00
4	☐			0.00
5	☐			0.00
6	☐			0.00
7	☐			0.00
8	☐			0.00
9	☐			0.00
10	☐			0.00

☐ Set factor on registering MC

☒ Base Shift

OK Cancel Apply Help

The diagram consists of three black arrows. The first arrow points from the 'None' radio button to the 'Disp.' column of the table below. The second arrow points from the 'Event' column of the table below to the 'Event' column of the table above. The third arrow points from the 'm/z' column of the table below to the 'm/z' column of the table above.

Method View - Peak Integration Parameters

View Edit

Integration Identification Quantitative Compound Performance Spectrum Library Custom QC Check

Mass/Compound Name List

Default
254.00(Ev1)

Integration

☒ Auto(Area) ☐ Auto(Height) ☐ Advanced

MaxPeak: 1

Width: 10 sec

Slope: 1000 /min

Drift: 0 /min

I. DBL: 1000 min

Min. Area/Height: 0 counts

Calculated by: ☒ Area ☐ Height

Program

Noise/Drift Calculation...

Smoothing

Method View - Quantitative Parameters

View Edit

Integration Identification Quantitative Compound Performance Spectrum Library Custom QC Check

Quantitative Method: External Standard

Unit: ppb

Calculated by: ☒ Area ☐ Height

Format of Concentration: ☒ ☐

Decimal Digits: 2 Significant Digits

Calibration Curve


of Calib. Levels: 3

Curve Fit Type: Linear

Zero: Not Forced


Weighting Method: None



Sample Information...




Double click để chọn
Ret. Time

Nhập nồng độ
của chuẩn.







Folder: F:\LCMS Data\23-12-2009

Filename	Program	Background	Background D	Sample Type	Analysis Type	Method File	Data File	Level#
BG001		✓	BG001.lcd	1:Standard(I)	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	Data001.lcd	1
Dat		✓	BG001.lcd	1:Standard	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	Data002.lcd	2
D		✓	BG001.lcd	1:Standard	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	Data003.lcd	3
C		✓	BG001.lcd	0:Unknown	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	sulfa sample001.lcd	0
sulfa sample001		✓	BG001.lcd	0:Unknown	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	sulfa sample002.lcd	0
sulfa sample002		✓	BG001.lcd	0:Unknown	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	sulfa sample003.lcd	0
sulfa sample003		✓	BG001.lcd	0:Unknown	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm	sulfa sample004.lcd	0
sulfa sample004		✓	BG001.lcd	0:Unknown	MIT MQT ILT STT	SULFAMETHOXAZOLE.lcm		
Sulfa std001								
Sulfa std002								
Sulfa std003								
Sulfa std004								
Sulfa std005								

